


- For use with 10-stage PMTs that fit standard 14-pin sockets
- Linear output available from anode or tenth dynode
- Focus control for optimum performance



The ORTEC Model 266 Photomultiplier Base provides voltage distribution to essentially all 10-stage photomultiplier tubes (PMTs) that fit its standard 14-pin tube socket. It provides capacitively-coupled linear output signals from the anode and the tenth dynode that can be used in either timing or linear pulse height analysis systems. This arrangement allows the use of either polarity of output without an inverting amplifier. High-quality signals are maintained when these outputs are fed through 50- $\Omega$  terminated cables or directly into a linear preamplifier such as the ORTEC Model 113 Scintillation Preamplifier. The Model 113 output is fed into one of the ORTEC main amplifiers, where it can be either active-filter- or delay-line-shaped for analysis. The focus control on the Model 266 allows optimum adjustments for the best performance of the particular PMT.

The Model 266 is compatible with most standard 10-stage PMTs that fit standard 14-pin sockets, including those listed in Table 1.

Other compatible tubes may be determined by comparison with those that are listed (also see Figs. 1 and 2). Note that all photomultiplier tube specifications are given by the PMT manufacturer.

**Table 1. Compatible Photomultiplier Tubes.**

| ADIT    | Burle (formerly RCA) | Hamamatsu     | Phillips | Electron |
|---------|----------------------|---------------|----------|----------|
| B51B01  | 4900                 | PM55          | XP2202   | 9266K    |
| L51B01  | 5819                 | R208          | XP2203B  | 9272K    |
| V51B01  | 6342A                | R550          | XP2412B  | 9250K    |
| B51D01  | 6655A                | R594          |          | 9256K    |
| B51C01  | S83006E              | R877          |          | 9305K    |
| B76B01  | S83013F              | R878          |          | 9265K    |
| V76B01  | S83019F              | R1507         |          | 9269K    |
| B76C01  | S83020F              | R1512         |          | 9273K    |
| B89B01  | S83021E              | R1513         |          | 9274K    |
| B89C01  | S83022F              | R1612         |          | 9306K    |
| B89D01  | S83025F              | R1791         |          | 9390K    |
| B133D01 |                      | R1836         |          | 9275K    |
| B133C01 |                      | R1847-07      |          |          |
| V133B01 |                      | R1848-07 7696 |          |          |

# 266

## Photomultiplier Base

### Specifications

#### PERFORMANCE

**BLEEDER RESISTANCE** 1.5 M $\Omega$  total, tapped to provide proportional bias steps to successive tube elements.

#### CONTROL

**FOCUS** The voltage to the focus electrode in the tube is available as an external trim adjustment using a screwdriver potentiometer.

#### INPUTS

**POS HV** SHV connector accepts positive bias voltage to 2.5 kV maximum.

**PMT SOCKET** TRW 3B14. Fits JEDEC B14-38 PMT pin base (see Fig. 2).

#### OUTPUTS

**ANODE** BNC connector provides negative linear output through  $Z_0 \sim 1.1$  M $\Omega$ , capacitively-coupled.

**DYNODE** BNC connector provides positive linear output from the tenth dynode through 1.1 M $\Omega$ , capacitively-coupled.

#### ELECTRICAL AND MECHANICAL

##### WEIGHT

**Net** 0.37 kg (0.81 lb).

**Shipping** 0.96 kg (2.12 lb).

**DIMENSIONS** 5.6 cm (2.2 in.) diam x 10.2 cm (4 in.) long.

used between the preamplifier output and the input to the main shaping amplifier. For timing from the anode use a C-25-12 Cable. The C-36-12 Cable is recommended for the high voltage connection.

### Ordering Information

To order the Model 266 Photomultiplier Base or related accessories, use the following model numbers and descriptions:

| Model   | Description   |
|---------|---|
| 266     | Photomultiplier Base  |
| C-24-1  | RG-62A/U 93- $\Omega$ Cable with two BNC male plugs, 1-ft length    |
| C-24-12 | Same as above, 12-ft length   |
| C-25-12 | RG-58A/U 50- $\Omega$ Cable with two BNC male plugs, 12-ft length   |
| C-36-12 | RG-59A/U 75- $\Omega$ Cable with two SHV female plugs, 12-ft length |
| C-27    | Terminator, 100 $\Omega$ , BNC male plug                            |




Fig. 2. JEDEC B14-38 PMT Pin Base, with Pin Assignments:

|          |                     |
|----------|---------------------|
| d1 - d10 | dynodes 1 to 10     |
| a        | anode               |
| i.c.     | internal connection |
| g        | grid                |
| k        | cathode             |




Fig. 1. Simplified Schematic Diagram of ORTEC Model 266 Photomultiplier Base.

Specifications subject to change  
043008

**ORTEC**<sup>®</sup>

[www.ortec-online.com](http://www.ortec-online.com)

Tel. (865) 482-4411 • Fax (865) 483-0396 • [ortec.info@ametek.com](mailto:ortec.info@ametek.com)  
801 South Illinois Ave., Oak Ridge, TN 37831-0895 U.S.A.  
For International Office Locations, Visit Our Website

**AMETEK**<sup>®</sup>  
ADVANCED MEASUREMENT  
TECHNOLOGY