
OPERATING INSTRUCTIONS

Version 3.7

polystat compatible control 1

Foreword: Please read carefully!	1
Product description	1
Safety Advice Regarding Overtemperature protection device	1
Safety Advice Regarding Controller Replacement	1
New functions from V3.04	2
G. Description of unit	3
G.1 Application	3
G.2 Preparation	4
G.3 Commissioning	4
G.4 Level alarm	5
G.5 Electrical thermostat within pump motor	5
G.6 Overtemperature protection device	5
G.7 Heat transfer fluid	5
G.8 Unit Diagram	6
G.9 Type recognition	6
G.9.1 Table of Device Types	6
G.9.2 Table of Controller Types	7
G.10 Operation of Refrigerated Baths K6, K12, K15, K20, K25	7
G.10.1 Service, maintenance	8
G.10.2 Transport	8
1. Features of the set-key	9
1.1 Reading the set-point value	9
1.2 Adjusting the set-point value	9
2 Description of the Menu Functions	10
2.1 Call up the Menu functions	10
2.2 LIMIt Setting set-point value limits	10
2.3 dEG. Setting the temperature-unit to Celsius or Fahrenheit	11
2.4 COMP. Set Automatic Compressor Control	11
2.5 HEAt Set Maximum Heating Capacity	11
2.6 rOM Read Software version	11
2.7 Init Set factory-set unit parameters	12
2.8 SHIFt Pt100 internal Offset (Shift) circuit calibration program	12
2.9 CAL Pt100 internal circuit calibration program	13
2.10 P_UP The automatic cutout	16
2.11 SAFE The backup memory for equipment configuration and calibration data	16
2.12 PUMP The Venting Function	18
2.13 Contr. Set parameters	19
2.14 tEMP. Start Thermoregulation	19
3 Alarm messages Polystat Compatible Control	20

Foreword: Please read carefully!

Product description

Polystat cc1 (61) is a compact, microprocessor control unit for Huber's Polystat Compatible Control circulation baths. The set-point limits, the compressor automatic controls and the °C or °F display can be set. The controller can be calibrated to comply with ISO9000ff. All important data is stored in a permanent memory. Data entry is controlled by the operator and supports alpha-numeric characters.

Safety Advice Regarding Overtemperature protection device

The overtemperature protection device is an independent safety device which must be set by the user. When setting this unit, allow for the flash-point of the heat transfer fluid and the flash-point of other media which are used in the application as well as the temperature limits for the application apparatus, i.e., the connecting hoses. Regularly (e.g., every 4 weeks) check the operation of the overtemperature protection device at a non-critical temperature. The user is responsible for ensuring that the setting of the overtemperature protection device is correct.

Polystat 61

The polystat 61 is available with or without overtemperature protection alternatively. You can recognize the version without overtemperature protection by the missing variable-pitch screw (position 4, fig. on page 5). The version without overtemperature protection is only allowed to be used with **non-flammable** liquids.

Safety Advice Regarding Controller Replacement

Electrical hazards must never be underestimated. For this reason, always switch off the circulator bath and disconnect the mains power supply before replacing the controller (21). Ensure that the mains power is not reconnected. Remove the fastening screws (22), take out the controller (21) forward. Ensure foreign bodies do not enter the control unit slot once the control unit has been removed from the circulator bath. Do not touch the control unit or the circulator bath contacts under any circumstances. Examine the replacement control unit contacts before installing the unit. The contact strip must not be dirty. When inserting the replacement control unit, ensure it is not tilted and do not use excessive force under any circumstances. Re-insert the fastening screws (22). The circulator bath can now be re-commissioned.

The data entry modes are divided into two levels.

1. Feature selection:

Features as set-point value adjustment, etc. are selected using the **set**-key. The display does **not** flash during these operations. The selected feature can be accessed using the set-key, thereby permitting entry to the 2nd programming level.


2. Data entry mode:

The adjuster is used to set a numerical value, e.g., the set-point value or to select a programme option, e.g., Auto, On or Off for the compressor automatic controls. The display flashes during this operation. The value or the option must then be confirmed using the set-key. The controller initially stores the data in a non-volatile memory. The control unit will annul the data entry and return to the default status (displaying the actual value) unless the data entry is confirmed within 4 minutes.

New functions from V3.04

From software version V3.04, a few additional functions have been incorporated in the operating software.

The thermostat is now fitted with a programmable automatic cutout. The user can determine the activation conditions himself. Only when the automatic cutout is switched on does the thermostat start the temperature control when the power is connected or after a power failure. If the automatic cutout is switched off, the user

has to start the temperature control when  is displayed by turning the rotary knob (see §2.10).

The definition of the calibration temperatures and the two-point calibration programme have been summarized under point CAL in the function menu (see §2.9). Additionally, you have the possibility to equalize the internal Pt100 sensor at a temperature. This Offset (Shift) calibration is described under point §2.8.

It is now possible to save the user-defined settings in the permanent memory so they may be reloaded if required. Access to this function is a complicated process in order to prevent it from being called up unintentionally (see §2.11).

Another innovation is the saving of the calibration settings. The user-saved calibration data or the factory calibration data may be loaded if required. Access to this function is a complicated process in order to prevent it from being called up unintentionally (see §2.11).

To improve external ventilation, a ventilation function has now been installed. This enables the pump to be operated in alternating mode. (see §2.12).

G. Description of unit

The Polystat is a temperature control unit with Plug & Play technology featuring microprocessor control and straightforward operation with an angular sensor and a digital display. High-pressure pump, coated components manufactured from stainless steel or plastic which is highly resistant to chemical attack. All models with 3 year guarantee for electronic components. Selection of 3 Models:

Polystat cc1, with level safety device (float switch) and adjustable overtemperature protection device for long-term unattended operation with flammable and non-flammable fluids (FL). Additional safety through maximum and minimum set-point value. Max. temperature=200°C*.

Polystat 61, as Polystat cc1, however alternatively with or without overtemperature protection. The version without overtemperature protection is only allowed to be used with non-flammable liquids (NFL). Max. temperature=100°C.

Polystat cc2 (62), as Polystat cc1 (61), however with program controller (5 steps), master control unit for external temperature control, 3 set-point temperatures and straightforward ramp feature. Max. temperature = 200°C*.

Polystat cc3 (63), as Polystat cc2 (62), however with the entire range of features from the Compatible Control circulation baths: Program controller (50 steps can be divided into 10 programmes), RS232, RS 485 and analogue interfaces (4...20mA) for two-way communication. Connection for serial printer. Max. temperature = 200°C*.

The **Polystat as an immersion temperature control unit** with a 2000 watt heater output which provides the features of an immersion circulation bath or a refrigerated circulation bath when combined with an immersion bath or a refrigerated bath. The pump permits temperature control of sealed equipment, e.g., photometers when used in conjunction with the pump adapter.

The **Polystat as a bath thermostat** (Heating thermostat). Maximum operating temperature see nameplate. For works at room temperature the cooling coil can be operated with water. With lower temperatures a separate cooling is needed.

G.1 Application

Temperature range and safety classification.

Immersion Temperature Control Units	Min. Temp.	Max. Temp.	Safety Classification	
Polystat CC 1 (61)	-30°C	200°C (100°C)	FL (NFL)	S
Polystat CC 2 (62)	-30°C	200°C	FL	S
Polystat CC 3 (63)	-30°C	200°C	FL	S
Heated Temperature Control Units				
Polystat 201	-30°C	200°C	FL	S
Polystat 202	-30°C	200°C	FL	S
CC 302	-20°C	300°C	FL	D
CC 303	-20°C	300°C	FL	D
Heated/ Refrigerated Temperature Control Units				
Polystat K6	-25°C	150°C	FL	S
Polystat K12	-20°C	120°C	FL	S
Polystat K15	-20°C	200°C	FL	S
Polystat K20	-30°C	200°C	FL	S
Polystat K25	-30°C	200°C	FL	S
Intelligent Chiller (-H), Unichiller UC (-H)	see nameplate	40°C (80°C)	NFL (FL)	D

S= Float switch, D= Flow indicator

* The temperature range can be extended up to 300 °C in conjunction with a CC302, CC303.

	Min. Temp.	Max. Temp.
Polycarbonate baths A5 up to A18	-30°C	100°C
Stainless steel baths insulated with PUR housing E8, E12, E20	-30°C	120°C
Stainless steel baths insulated with stainless steel housing E15, E25	-30°C	200°C
Refrigerated baths:		
Polystat K6	-25°C	(80°C) / 150°C
Polystat K12	-20°C	(50°C) / 120°C
Polystat K15	-20°C	(50°C) / 200°C
Polystat K20	-30°C	(50°C) / 200°C
Polystat K25	-30°C	(50°C) / 200°C

(50°C) = Maximum working temperature with refrigeration unit switched on.

For polycarbonate and stainless steel vessels: A water-filled refrigeration coil is required to cool the fluid below its inherent temperature. An external mechanical refrigeration system must be fitted where even lower fluid temperatures are required. There are basically two options for external refrigeration: 1st an immersion refrigeration unit with refrigeration unit sensor located directly temperature-controlled bath or 2nd a pipework refrigeration unit located at the pump adapter.

G.2 Preparation

Note the following information before commissioning

- Supply voltage and fusing - see name plate.
- Maximum working temperature - see tables in § G.1. In addition, observe the working temperatures of the heat transfer fluid (viscosity and flash-point).
- A immersion temperature control unit is fastened to the temperature control vessel using the clamping bolt (20). This can be adjusted according to the application. To do this, remove the three fastening bolts for the bolted clamp and fasten the bolted clamp in the desired position.
- Max. bath filling level = 15mm below the top edge of the bath
Please note: Allow for a volumetric expansion of approx. 15% when heating the heat transfer oil from 20 °C up to 200 °C.
- Connect the external equipment to the pump adapter (accessory at present) or seal the equipment using the closure bolts. An SW19 open-ended spanner is required for the union nut - counterhold using an SW 17 open-ended spanner. Always counterhold when connecting or disconnecting the equipment!
- Connect the coolant water hose if the fluid temperature has to be regulated to a temperature below its inherent temperature. An immersion refrigeration unit can be used to achieve colder temperatures, e.g., HUBER TC40. An equivalent refrigeration effect can be achieved using a pipework refrigeration unit, e.g., HUBER DC30. This is connected to the pump adapter (15 + 16). The return line from the equipment is routed via the pipework refrigeration unit if external temperature control is used.

G.3 Commissioning

Check the level of the bath and all hose connections. Connect the mains power supply. Activate the temperature control unit using the mains power switch (2). *HELLO* will appear on the display (5) followed by the temperature control unit identification. Consult the following sections for information regarding operation of the controller.

G.4 Level alarm

The float switch (9) breaks the circuit before the fluid level falls below a level which is critical for the heater (12). The message *Float* will appear on the display (5) and the temperature control unit will deactivate. The *Float* message will be retained, even if the cause of the message is rectified (refilling with heat transfer liquid). For cancelling the message, deactivate and then reactivate the temperature control unit.
Max. bath filling level = 15 mm below the top edge of the bath. Please note the interlock temperature of the volume.

G.5 Electrical thermostat within pump motor

If the pump motor overheats, the motor's internal electrical thermostat will trip out. This may result from excessive heat transfer fluid viscosity. A message will appear on the display.

With units with float switch the message *pump* appears: the thermostat switches off.

With units with flow indicator the message *float* appears: because of failure of the pump motor the pressure switch responded.

The message will be retained even once the pump motor has cooled down. This can process take some time, depending on the bath temperature. For cancelling the message, deactivate and then reactivate the temperature control unit.

G.6 Overtemperature protection device

For devices with security classification FL.

If the bath temperature exceeds the value set for the overtemperature protection device, the *Temp* message will appear on the display and the temperature control unit will deactivate. The *Temp* message will be retained even once the bath has cooled down. This procedure can take some time, depending on the bath temperature. For cancelling the message, deactivate and then reactivate the temperature control unit.

G.7 Heat transfer fluid

Ensure the heat transfer liquid used for temperature control is either water or one of the fluids specified in our price list.

If using water as a heat transfer fluid, ensure that the freezing point of water is accounted for.

Compare the minimum and maximum permissible temperatures for the heat transfer fluid against your desired range of working temperatures. The maximum working temperature should remain at least 5K below the flash-point of the heat transfer fluid. Please note that the viscosity of the heat transfer fluid increases as its temperature is reduced. With the minimum operating temperature the viscosity of the thermo fluid must not exceed 50 mm²/s. The freezing point should stay 10-20K below the operating temperature. Likewise, check the controller settings (set-point limits §2.2).

Use of the following fluids is not permitted:


- demineralised or distilled water
- mineral-enriched water or sea water
- CaCl₂ salt solution
- thermoregulation fluids containing ether, ester and amine. These admixtures are in some sorts of glycol as these heat transfer fluids can damage the 1.4301 stainless steel construction material.

Please close the pump sockets with screw plugs, resp. with external application connect the tubes. The tube material must be suited to the requested operating temperature range. Verify pump sockets to density.

G.8 Unit Diagram

- 1 Hood
- 2 Mains power switch
- 3 Mains power connection
- 4 Overtemperature protection device
- 5 Display
- 7 Rotary knob
- 8a Power LED
- 8b Heating LED
- 8c External LED
- 8d Program LED
- 9 Float (level)
- 10 Controller sensor

- 11 Stirrer
- 12 Heater
- 13 High-pressure pump
- 14 Fastening for 15-16
- 15 Pump adapter pressure connection 6 (accessories)
- 16 Pump adapter return connection (accessories)
- 18 Cooling coil pipework connection (accessories)
- 19 Fastening for 18
- 20 Bolted clamp
- 21 Controller
- 22 Fastening screw
- 23 6pole connector

G.9 Type recognition

Device type recognition and sensor type recognition first featured on the Polystat. As a result, this permitted the use of parameters which are specific to the device and the controller. The currently defined devices and controllers are described in the following table. The „TYPError„ fault message is displayed if a controller with a unrecognised device type. The device table can be revised by updating the software. The *HELLO* display indicates that a device has been recognised.

G.9.1 Table of Device Types

Device	Display		Note
Polystat NFL	HELLO	P_100	Without overtemperature protection device, -30°C... + 100°C
Polystat FL	HELLO	P_200	With overtemperature protection device, -30°C... + 200°C
Polystat 300°C	HELLO	P_300	With overtemperature protection device, -30°C... + 300°C
Intelligent Chiller Unichiller UC	HELLO	IC1...4	Without overtemperature protection device, -25°C... + 40°C
IC...-H Unichiller UC...-H	HELLO	IC1-H...4-H	With overheating protection device, -25°C... + 80°C

G.9.2 Table of Controller Types

Controller	Display	Note
Polystat CC1 (61)	Contr. 1/ Contr. 2	Set-point value limit / extended temperature range
Polystat CC2 (62)	Contr. 2	As Poly CC1 (61) plus external circuit temperature controller, 3 user temperatures, small programme controller, straightforward ramp procedure
Polystat CC3 (63)	Contr. 3	As Poly CC2 (62) but 10 user temperatures, larger program controller plus RS232 interface, RS485 interface, 4..20mA analogue interface

G.10 Operation of Refrigerated Baths K6, K12, K15, K20, K25

Water baths without approaching possibility of regulator : With these water bathes no 6-pole connector at the rear of the bath.

Connect the mains power supply using the heater unit plug on the back of the housing. One of two temperature control unit operation modes can be selected.

- Cooling is not required or the operating temperature is higher than the maximum working temperature of the cooling machine:

Move the black rocker switch to the 0 position (located above the voltage supply unit on the back of the housing). Activate the temperature control unit using the mains power switch (2). The circulation pump and the temperature control unit are activated.

- Cooling is required:

The commissioning procedure is performed as described above, however the refrigeration unit is also activated using the black rocker switch on the back of the housing. Please note that the refrigeration unit can only be operated up to the maximum working temperature (see table §G1).

The refrigeration unit must be deactivated for bath temperatures above the maximum value using the black rocker switch on the back of the housing.

Please ensure that the bath temperature and the set-point value are indicated on the LED display once the Polystat is activated.

Water baths with approaching possibility of regulator : With these water bathes a 6-pole connector is at the rear of the bath.

Connect the data to the regulator. Use the delivered connection cable. Set up the power supply with the high temperature plug at the rear of the unit. Only change the rocker switch at the rear of the water bath to position 1 if the bath is working without data control. Then continue as described above. If not the rocker switch has to be in position 0.

You have the possibility to choose between two ways of operation of the thermostat.

- You do not need a cooling capacity or the operating temperature is higher than the maximum working temperature of the cooling machine.

Switch on the thermostat with the main switch (2). The circulation pump and the tempering have started. Switch off the cooling machine with the compressor automation (see §3.11)

- You need cooling capacity and the operating temperature is lower than the maximum working temperature of the cooling machine:

For installation see above, however the automation of compressing can be switched to AUTO or ON (see §3.11) Please see that the operation of the cooling machine is possible only up to a maximum working temperature (see table §G1). With bath temperatures above the maximum working temperature the cooling machine has to be switched off through the compressor automation (see §3.11).

In case of simultaneous operation of cooling bath (compressor) and polystat (heating), the heating capacity has to be limited to a maximum of 75% in order not to exceed the maximum current absorption of 10A.

Therefore, the function *Heat* has been introduced in the function menu. With this function, the maximum heating capacity can be set on 25%, 50%, 75% or 100% (see §3.12).

If the heating capacity is not limited, the fuse integrated in the cooling bath's bush is activated. The bush is at the cooling bath's back.

G.10.1 Service, maintenance

Ensure there is an approx. 10 mm gap on both sides of the unit.

Clean condenser regularly. Remove the cover plate by releasing both of the fastening bolts and clean the condenser with compressed air.

G.10.2 Transport

The thermostats do have to be secured for transport, however the contents of the bath should be drained completely. The temperature control unit should not be tilted or placed on its side when transported.

Ensure the unit is well packaged and mark the packaging with the correct vertical transportation orientation.

1. Features of the set-key

The set key is for choosing the setpoint or menu functions and, in general, for confirming an input.

Safety first!

The set-point value limit should be adjusted when commissioning the unit and each time the temperature regulation liquid is replaced.

The set-point value limit is software-controlled and separate from the high-temperature and high-level safety devices whose activation does is not controlled by the control system electronics.


The setpoint limit is no overtemperature protection.

Limit the upper set-point value to 5K below the flash-point of the temperature regulation liquid (e.g., 25°C when ethanol is used).



Limit the lower set-point value according to the viscosity or freezing point of the temperature regulation liquid (e.g. 5°C with water).

See section §2.2 for the relevant procedure.

1.1 Reading the set-point value

   	Press the set-key once, then release The adjusted set-point value is displayed in the units selected (°C/°F) and flashes. e.g.: 25.0° Centigrade or or 77.0 ° Fahrenheit After approx. 4 minutes the actual value is displayed.
---	---

1.2 Adjusting the set-point value

   	Press the set-key once. The previously saved set-point is displayed flashing. Use the rotary knob to raise or lower the set-point value, as required. Press the set-key once to confirm the set-point value. If the set-point value is not confirmed within 4 minutes, the data entry mode will be terminated and the adjusted set-point value will not be saved.
 	

2 Description of the Menu Functions

2.1 Call up the Menu functions

The function menu can only be called up when starting the thermostat. In this menu all other functions of the thermostat can be set.

   5 seconds      	Press the set key and switch on the thermostat at the mains switch. Wait is displayed. Release the set key. After approx. 5 Sec. the first function of the function menu <i>Limit</i> is displayed. You can leaf through the different functions of the function menu using the rotary knob. By pressing the set key, the chosen function is called up. You can exit the function menu when <i>Temp</i> is displayed. The thermoregulation starts by choosing the function with the help of set key. The bath temperature is displayed.
--	--



2.2 LIMIt Setting set-point value limits

Set the maximum set-point value to at least 5 K below the flash-point or boiling point (water). Set the minimum set-point value to at least 3 K above the freezing point of the temperature regulation liquid and taking account of its viscosity. Generally, increased viscosity results in reduced heat transfer rate and thereby reduced system efficiency.

The upper and lower limits of the set-point value should be redefined each time the temperature regulation liquid is replaced in order to ensure operation within a suitable temperature range and also accounting for the flash-point and viscosity. If flammable liquids are used (e.g., alcohol), the upper limit must be 5K below the flash-point.

             	In case of need call up the function menu (see 2.1). Turn the rotary knob until <i>Limit</i> is displayed. Once the set key is pressed <i>L-Lim</i> will be initially displayed. The lower set-point value limit can now be adjusted using the rotary knob. Example: The lower limit is 2°C Confirm the temperature setting using the set-key. <i>H-Lim</i> is displayed. The upper set-point value limit can now be adjusted using the arrow-keys. Example: The upper set-point value limit is 95°C. Once the settings are confirmed using the set-key. The function menu is still activated.
---	--

2.3 dEG. Setting the temperature-unit to Celsius or Fahrenheit

   dEG	In case of need call up the function menu (see 2.1). Turn the rotary knob until <i>Deg.</i> is displayed. Select the function with the Set key.
  dEG°C  dEG°F	Switch over between Deg°C and Deg°F using the rotary knob. Confirm the temperature measurement unit setting using the set-key. The function menu is still activated.

2.4 COMP. Set Automatic Compressor Control

   COMP	In case of need call up the function menu (see 2.1). Turn the rotary knob until <i>Comp.</i> is displayed. Confirm with set key.
  AUTO  ON  OFF	Choose the required compressor mode using the rotary knob. Confirm with set key. The function menu is still activated.

2.5 HEAt Set Maximum Heating Capacity




  HEAt	In case of need call up the function menu (see 2.1). Turn the rotary knob until the message <i>Heat</i> appears. Confirm with set key.
  100  75  50  25 	Set the required maximum heating capacity in % using the rotary knob. See §G.10 „Operation of Refrigerated Baths„. Confirm with set key. The function menu is still activated.

2.6 rOM Read Software version

  rOM	In case of need call up the function menu (see 2.1). Turn the rotary knob until <i>Rom</i> is displayed. Confirm with set key.
 0302 P_200 rEOL 1 2	Successively, the following is displayed: - the software version e.g. 03.02, - the type of device e.g. P_200 - the message <i>Regl.</i> for controller and finally the controller type, e.g. 1 The function menu is still activated.

2.7 Init Set factory-set unit parameters

In case of factory-set parameters, the control parameters and the set-point limits are cleared.





        	In case of need call up the function menu (see 2.1). Turn the rotary knob until <i>Init</i> is displayed. Confirm with set key. Yes is displayed. In order to set the factory-set parameters, confirm with set key. If you do not wish to change the parameters, set the display <i>No</i> with the rotary knob and confirm with set key. The function menu is still activated.
---	--

2.8 SHIFt Pt100 internal Offset (Shift) circuit calibration program

This program only starts when the internal circuit control mode is activated (see §3.2).

The internal circuit Pt100 sensor can be calibrated using the keyboard thus permitting calibration of the microprocessor control unit and the circulator bath for conformance with ISO 9000 ff. or for quality assurance requirements. The internal circuit sensor is calibrated at any point within the temperature range of the device. The Shift calibration can be adjusted a maximum of +/-5K.

         	In case of need call up the function menu (see 2.1). Turn the rotary knob until <i>Shift</i> is displayed. Call up the Shift program using the set key. The current Shift is displayed, e.g., -1.55K. The Shift is added to the temperature value. E.i. the indicated temperature rises in case of a positive shift and goes down in case of a negative shift. Use the rotary knob keys to adjust the Shift. Confirm the entry using the set-key. Finally, the shift setting is displayed for 2 seconds. The function menu is still activated.
--	--

2.9 CAL Pt100 internal circuit calibration program

The internal circuit Pt100 sensor can be calibrated using the keyboard thus permitting calibration of the microprocessor control unit, the bath and circulation circulator bath - thereby complying using ISO 9000 ff. or quality assurance requirements. The settings are saved once the calibration procedure is completed.

The control unit for the circulator bath is calibrated ex works. If the control unit is replaced, the system accuracy should be checked.

The calibration temperatures must be defined in the function menu (see §2.8). The second calibration temperature must always be higher than the first calibration temperature. The difference between the calibration temperatures must be greater than 10K and less than 250K.

Before performing the calibration procedure, the set-point limit values must be defined to ensure that they do not conflict with the calibration temperature. (see §2.2)

The display for the type of temperature regulation liquid depends on the calibration temperatures. At a calibration temperature of <20°C ALCOH. is displayed, between 20°C and 94°C H-2-O for water is displayed and at a calibration temperature of >=95°C OIL is displayed. These are only suggestions. The correct temperature regulation liquid must always be selected by the user.

The following describes the calibration of an internal circuit sensor at 0°C with alcohol and at 80°C with water:





  CAL	Fill the bath with alcohol and switch on the thermostat. Call up the function menu (see 2.1). Turn the rotary knob until <i>Cal</i> is displayed and confirm with set key.
Set  CAL_1	<i>Cal_1</i> appears on the display.
Set  CC 0	To input the first calibration temperature turn the rotary knob.
Set  OTHER	The rotary knob may be used to select from a table with predefined calibration temperatures, e.g. CC0 for 0°C. With the display <i>Other</i>, there is an option of entering any calibration temperature.
Set  CC 25	Confirm the table value or the option <i>Other</i> with the Set key. Only with the option <i>Other</i>:
Set  CAL_2	Use the rotary knob to set any calibration temperature and confirm with the Set key.
Set  CC 80	To enter the second calibration temperature turn the rotary knob when the display shows <i>CAL_2</i>.
Set  OTHER	The rotary knob may be used to select from a table with predefined calibration temperatures, e.g. CC80 for 80°C.
Set  CC 60	With the display <i>Other</i>, there is an option of entering any calibration temperature.
Set  START	Confirm the table value or the option <i>Other</i> with the Set key. Only with the option <i>Other</i>:
Set  no	Use the rotary knob to set any calibration temperature and confirm with the Set key.
Set  YES	<i>Start</i> appears on the display. To start the calibration program, select the option <i>Yes</i> with the rotary knob and confirm with set key. With the option <i>no</i>, the function is left.
 ALCOH	After the option <i>Yes</i>, the display <i>Alcoh</i> appears.
 CC 0 20.00	Acknowledge with the rotary knob. The program <i>CC 0</i> starts automatically, the controller specifies a set value of 0°C. The display <i>CC 0</i> appears for about 1s and alternates for 2s with the actual temperature eg: 20.00(°C).
 ENTER 0.00 0.00	If 0°C is measured and the bath temperature is stable for a certain time, the message <i>Enter</i> appears. Check the stability of the bath temperature on the reference thermometer. If this is stable, set the temperature shown on the reference thermometer with the rotary knob. The displayed value may be corrected by $\pm 5^\circ\text{K}$. Pressing the Set key confirms the entered temperature.
Set H-2-O	The display <i>H-2-O</i> for water appears.

2.10 P_UP The automatic cutout

The automatic cutout determines the behaviour of the thermostats when the power is switched on or after a power outage. If the automatic cutout is activated, the thermostat starts the temperature control when it is switched on at the mains switch. The thermostat also starts the temperature control again automatically after a power outage.

For safety reasons, it may be necessary for the thermostat to remain switched off after a power outage. In such a case, the automatic cutout must be switched off.

		In order to set the automatic cutout, press the set key when <i>P_UP</i> is displayed.
		Choose <i>Yes</i> for automatic cutout On or the option <i>No</i> for automatic cutout OFF with the help of the rotary knob.
	 	Confirm with set key.
	 	When the automatic cutout is switched off, the thermostat does not automatically start the temperature control. After the starting phase when the display shows <i>HELLO</i> and the type display, the display alternates between the current actual value and the reference <i>turn</i> . If the rotary knob is now turned, the thermostat starts the temperature control. The reference <i>turn</i> no longer appears on the display.

2.11 SAFE The backup memory for equipment configuration and calibration data

The access to this function is a complicated process to prevent it from being called up unintentionally. It is now possible to save the user-defined settings in the permanent memory so that it may be reloaded if required. This is displayed by the menu points *L.USER* and *S.USER*. All settings and programming operations are saved and loaded.

Another innovation is the saving or reloading of the calibration settings. The user-saved calibration data are saved under the menu point *S_CAL* and reloaded under the menu point *L_CAL*. The factory calibration data may be loaded under the menu point *ReCAL* if required. Access to the load function is protected with a *Yes/No* question to prevent it from being called up unintentionally.

The calibration data should only be reloaded if the user is aware of the content of the backup copy. If false data are backed up, loading this data may degrade the accuracy.

	SAFE	In order to access the main memory for the unit configuration and the calibration data, press the set key when <i>Safe</i> is displayed. The reference <i>Code</i> appears on the display.
	Code	
	108	Use the rotary knob to sent the password <i>108</i> and confirm with the Set key.
	LUSER	Use the rotary knob to set the desired backup option. <i>L.User</i> stands for load user data (equipment configuration),
	SUSER	<i>S.User</i> stands for backup user data (equipment configuration),
	L_CAL	<i>L_CAL</i> stands for load user calibration data,
	S_CAL	<i>S_CAL</i> stands for backup user calibration data,
	ReCAL	<i>ReCAL</i> stands for load factory calibration,
	Quit	When the display shows <i>Quit</i> , the menu may be left without changing the data.
	POP U	If <i>L.User</i> is confirmed with the Set key, <i>POP U</i> appears on the display. The last equipment configuration saved is loaded from the backup memory.
	PUSHU	If <i>S.User</i> is confirmed with the Set key, <i>PushU</i> appears on the display. The current equipment configuration is loaded into the backup memory.
	SURE	If <i>L_CAL</i> is confirmed with the Set key, <i>Sure</i> appears on the display.
	no	Use the rotary knob to set <i>Yes</i> and confirm with the Set key.
	YES	
	POP C	<i>POP C</i> appears on the display, the last saved calibration data are loaded from the backup memory into the main memory.
	PUSHC	If <i>S_CAL</i> is confirmed with the Set key, <i>PushC</i> appears on the display. The current calibration data are loaded into the backup memory.
	SURE	If <i>ReCAL</i> is confirmed with the Set key, <i>Sure</i> appears on the display.
	no	Use the rotary knob to set <i>Yes</i> and confirm with the Set key.
	YES	
	ReCAL	<i>ReCAL</i> appears on the display, the factory calibration data are loaded from the backup memory into the main memory.

2.12 PUMP The Venting Function

In case of an external operation, an alternating operation of the pump can be very helpful when venting. The On- and OFF-time of the pump can be defined successively.

      	In order to activate the venting function press the set key when <i>Pump</i> is displayed. <i>Pu On</i> is displayed for the On-time. Set the required On-time of the pump with the help of the rotary knob and confirm with set key. The On-time range can be set between 3...90 seconds. <i>PuOff</i> is displayed for the Off-time. Set the required Off-time of the pump with the help of the rotary knob and confirm with set key. The Off-time range can be set between 3...90 seconds. During the On-time the actual value and <i>PU xx</i> are displayed alternatingly, xx refers to the remaining On-time. The power LED is off. During the Off-time the actual value and <i>Offxx</i> are displayed alternatingly, xx refers to the remaining Off-time. Terminate the venting function with the set key the thermoregulation starts and the power LED is on.
---	--

2.13 Contr. Set parameters

With the help of the P-I-coefficients you can influence the action of the thermostat.

       	In order to set the parameters, press the Set key when <i>Contr.</i> is displayed. Afterwards, the Shift is displayed (see §2.8). Select a pre-defined parameter set (Pid 0 ... Pid 8) with the help of the rotary knob. If the parameter setting PID 9 is selected, the freely programmable parameters are taken over. If <i>MEM</i> is selected, the parameters for parameter set PID 9 (proportional band and integral time parameters) can be set. <i>ProP</i> for proportional coefficient will be displayed. Now set the coefficient with the help of the rotary knob. <i>InteG</i> for integral coefficient is displayed. Now set the coefficient with the help of the rotary knob. Set with the Set key.
--	--

Pre-defined parameter sets are:

Pid.	0	1	2	3	4	5	6	7	8	9
Prop.	5000	3000	10000	5000	200	200	500	3000	5000	50...30000
Integr.	1000	1000	1000	500	200	0	0	0	0	0...30000

2.14 tEMP. Start Thermoregulation

  	The function menu is activated (see §2.1). In order to leave the function menu, turn the rotary knob until the message <i>Temp</i> is displayed. Confirm with set key. Thermoregulation starts, the type of unit is displayed.
---	---

3 Alarm messages Polystat Compatible Control

FLoAR	Level in bath too low.
EENP	Mech. overtemperature protection device has tripped
PUNP	Winding thermostat of the pump motor has responded.
PtInT	Internal circuit Pt100 sensor faulty or not connected
PtOut	External circuit Pt100 sensor faulty or not connected
PtCon	Condensation sensor faulty or not connected
TYPEr	The control unit is installed in a module which is not recognised by the control unit.
H-Err	Error in controller hardware. Inform Service ☎.
EEPRr	EEPROM error, notify Service ☎.
UnErr	Unplug the unit, notify Service ☎. Exchange Polystat Compatible Control
Turn	Automatic cutout not activated. Turn the rotary knob to start it.
Cond	In addition, the <i>Cond</i> message is displayed as soon as the condensation temperature exceeds 50°C. The cooling capacity is reduced. The alarm does not lead to deactivation of the unit.

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New Immersion Circulators with Plug&Play Technology - microprocessor controlled, easy handling with encoder and digital display. Powerful force pump, moistened parts in stainless steel or high-resistant plastics. All models with 3 years warranty for electronic. 3 models available:

Polystat cc1, with level protection (float switch) and adjustable overtemperature protection for continuous operation without personal assistance for using flammable and non-flammable liquids (FL). Maximum and minimum setpoint for additional safety. Max. temperature 200°C.

Polystat cc2 (62), similar to Polystat cc1, but with programmer (5 steps), temperature sequence controller for external thermoregulation, 3 fixtemperatures and easy ramping function. Max. temperature 200°C.

Polystat cc3 (63), similar to Polystat cc2, but with all functions of the Compatible Control Thermostats: Programmer (50 steps, divisible into 10 programmes) and interfaces RS232, RS485 and analog (4...20mA) for bidirectional communication. Serial printer output. Max. temperature 200°C.

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Accessoires* and periphery: Cooling coil for water cooling, pump adapter for external thermoregulation, adapter nom. dia 8 mm/ 12 mm, dummy plugs and sleeve nuts thread M16x1, micro boltings, connection tubes, connection cable*. Waterbathes in polycarbonat and stainless steel and units for cells and glasses.

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New Circulators with Plug&Play Technology - microprocessor controlled, easy handling with encoder and digital display. Powerful force pump, moisted parts in stainless steel or high-resistant plastics. All models with 3 years warranty for electronic. Cooling coil for (tap) water. Max. temperature 200°C. 3 models available:

Polystat 201-1, with level protection (float switch) and adjustable overtemperature protection for continuous operation without personal assistance for using flammable and non-flammable liquids (FL). Maximum and minimum setpoint for additional safety..

Polystat 201-2, similar to Polystat 201-1, but with programmer (5 steps), temperature sequence controller for external thermoregulation, 3 fixtemperatures and easy ramping function.

Polystat 201-3, similar to Polystat 201-2, but with all functions of the Compatible Control Thermostats: Programmer (50 steps, divisible into 10 programmes) interfaces RS232, RS485 and analog (4...20mA) for bidirectional communication (FL). Serial printer output.

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Accessoires* and periphery: adapter nom. dia 8 mm*/ 12 mm*, dummy plugs* and sleeve nuts thread M16x1*, connection tubes, connection cable.

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New Circulators with Plug&Play Technology - microprocessor controlled, easy handling with encoder and digital display. Powerful force pump, moisted parts in stainless steel or high-resistant plastics. All models with 3 years warranty for electronic. Cooling coil for (tap) water. Max. temperature 200°C. 3 models available:

Polystat 202-1, with level protection (float switch) and adjustable overtemperature protection for continuous operation without personal assistance for using flammable and non-flammable liquids (FL). Maximum and minimum setpoint for additional safety.

Polystat 202-2, similar to Polystat 202-1, but with programmer (5 steps), temperature sequence controller for external thermoregulation, 3 fixtemperatures and easy ramping function.

Polystat 202-3, similar to Polystat 202-2, but with all functions of the Compatible Control Thermostats: Programmer (50 steps, divisible into 10 programmes), interfaces RS232, RS485 and analog (4...20mA) for bidirectional communication (FL). Serial printer output.

Technische Daten	Technical Data	polystat 202-1-S3	202-2-S3	202-3-S3
Arbeitstemperaturbereich	Operating temperature range		70...200°C	
mit Wasserkühlung	with water cooling		20...200°C	
mit Kühlgerät	with refrigeration chiller		-30...200°C	
Temperaturkonstanz bei 70°C	Temperature stability at 70°C		0,05 K (DIN 58966)	
Temperatureinstellung	Temperature adjustment		digital	
Temperaturanzeige	Temperature indication		digital	
Absolutgenauigkeit	absolute accuracy	kalibrierbar / setup for calibration		
Temperaturfühler	Temperature sensor		Pt100	
externer Temperaturfühler	external Temperature sensor	---	Pt100	Pt100
externer Programmeingang	external program input	---	---	4-20mA
Schreiber Ausgang	Recorder output	---	---	4-20mA
Sicherheitsklasse	Safety classification		FL	
Heizleistung	Heating capacity		1,0 kW	
Druckpumpe Anschluß 12mm	Force pump (adapter nom 12 mm)		10 l/min max.	
Anschluß 8 mm	(adapter nom 8 mm)		7 l/min max.	
Förderhöhe (Druck)	Pressure		0,2 bar max.	
Pumpenanschluss	Pump connection		M16x1	
Pumpenanschluss für Schlauch	Pump connections for hose		NW8/ NW12	
Füllvolumen	Bath capacity		4,5 lit.	
mit Verdrängereinsatz	with displacement rack		3,1 lit.	
Badöffnung B x T / Tiefe	Bath opening w x l / depth		105 x 90 mm / 150 mm	
Abmessungen B x T x H	Overall dimension w x l x h		150 x 285 x 350 mm	
Arbeitshöhe Bad	Heith of bath opening		190 mm	
Gewicht, netto	Net weight		8,3 kg	
Netzanschluss	Power supply requirement		110-120 V ~ 50/60 Hz	
Leistungsaufnahme/Absicherung	Power input / fuse		2100 Watt 10 A	
Bestell-Nr.	Order-No.	686.0001-S3	686.0002-S3	686.0003-S3
ab Fert.Nr.	from serial no.			V1.0

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Polystat bath thermostat consisting of immersion circulators polystat cc, transparent baths and bath bridge. Baths made of polycarbonate with temperature range up to 100°C. Bath bridge for A11 and A18 with bore hole for cooling probe (e.g. for HUBER immersion cooler TC40-TC100E).

name:

example: A5-2 transparent bath with 5 l bath capacity, in combination with polystat cc2
A11-3 transparent bath with 11 l bath capacity, in combination with polystat cc3

Technische Daten	Technical Data	A5-1	A5-2	A5-3
Füllvolumen.	Bath capacity		5 lit.	
Badöffnung B x T / Badtiefe	Bath opening w x l / bath depth	120 x 110 / 150 mm		
Abmessungen B x T x H	Overall dimension w x l x h	132 x 275 x 320 mm		
Gewicht, netto	Net weight		4 kg	
Testglaseinsätze Typ A bis H	tube racks typ A to H		1 St.	
Bestell-Nr.	Order-No.	688.0011-S3	688.0012-S3	688.0013-S3

Technische Daten	Technical Data	A6-1	A6-2	A6-3
Füllvolumen.	Bath capacity		6 lit.	
Badöffnung B x T / Badtiefe	Bath opening w x l / bath depth	120 x 210 / 150 mm		
Abmessungen B x T x H	Overall dimension w x l x h	132 x 375 x 320 mm		
Gewicht, netto	Net weight		4 kg	
Testglaseinsätze Typ A bis H	tube racks typ A to H		2 St.	
Bestell-Nr.	Order-No.	688.0014-S3	688.0015-S3	688.0016-S3

Technische Daten	Technical Data	A7-1	A7-2	A7-3
Füllvolumen	Bath capacity		7 lit.	
Badöffnung B x T / Badtiefe	Bath opening w x l / bath depth	120 x 310 / 150 mm		
Abmessungen B x T x H	Overall dimension w x l x h	132 x 475 x 320 mm		
Gewicht, netto	Net weight		4,5 kg	
Testglaseinsätze Typ A bis H	tube racks typ A to H		3 St.	
Bestell-Nr.	Order-No.	688.0017-S3	688.0018-S3	688.0019-S3

Technische Daten	Technical Data	A11-1	A11-2	A11-3
Füllvolumen.	Bath capacity		11 lit.	
Badöffnung B x T / Badtiefe	Bath opening w x l / bath depth	292 x 165 / 150 mm		
Abmessungen B x T x H	Overall dimension w x l x h	308 x 331 x 325 mm		
Gewicht, netto	Net weight		6,5 kg	
Testglaseinsätze Typ 1 oder 2	tube racks typ 1 or 2		2 St.	
Bestell-Nr.	Order-No.	688.0020-S3	688.0021-S3	688.0022-S3

Technische Daten	Technical Data	A18-1	A18-2	A18-3
Füllvolumen	Bath capacity		18 lit.	
Badöffnung B x T / Badtiefe	Bath opening w x l / bath depth	292 x 325 / 150 mm		
Abmessungen B x T x H	Overall dimension w x l x h	308 x 491 x 325 mm		
Gewicht, netto	Net weight		7,5 kg	
Testglaseinsätze Typ 1 oder 2	tube racks typ 1 or 2		4 St.	
Bestell-Nr.	Order-No.	688.0023-S3	688.0024-S3	688.0025-S3

Zubehör	Accessoires:	A5	A6	A7	A11	A18
Entleerungshahn	drain cock	6026	6026	6026	6026	6026
Badbrücke	bath bridge	6309	6309	6309	6310	6310

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Polystat bath thermostat consisting of immersion circulators polystat cc, insulated stainless steel bath with PUR-housing (B8,B12,B20), resp. with housing made of stainless steel (B15, B25) and bath bridge. Bath bridge with bore hole for cooling probe (e.g. for HUBER immersion cooler TC40-TC100E), and

name:

example: B8-2 stainless steel bath with 8 l bath capacity, in combination with polystat cc2
B25-3 stainless steel bath with 25 l bath capacity, in combination with polystat cc3

Technische Daten	Technical Data	B8-1	B8-2	B8-3
Arbeitstemperaturbereich	Operating temperature range		max 120°C	
Füllvolumen.	Bath capacity		8,5 lit.	
Badöffnung B x T / Badtiefe	Bath opening w x l / bath depth	235 x 168 / 150 mm		
Abmessungen B x T x H	Overall dimension w x l x h	290 x 350 x 375 mm		
Gewicht, netto	Net weight	11 kg		
Bestell-Nr.	Order-No.	688.0026-S3	688.0027-S3	688.0028-S3

Technische Daten	Technical Data	B12-1	B12-2	B12-3
Arbeitstemperaturbereich	Operating temperature range		max 120°C	
Füllvolumen	Bath capacity		12 lit.	
Badöffnung B x T / Badtiefe	Bath opening w x l / bath depth	300 x 193 / 150 mm		
Abmessungen B x T x H	Overall dimension w x l x h	355 x 380 x 375 mm		
Gewicht, netto	Net weight	15 kg		
Testglaseinsätze Typ 1 bis 4	tube racks typ 1 to 4	2 St.		
Bestell-Nr.	Order-No.	688.0029-S3	688.0030-S3	688.0031-S3

Technische Daten	Technical Data	B15-1	B15-2	B15-3
Arbeitstemperaturbereich max	Operating temperature range max		200°C	
Füllvolumen.	Bath capacity		15 lit.	
Badöffnung B x T / Badtiefe	Bath opening w x l / bath depth	300 x 193 / 200 mm		
Abmessungen B x T x H	Overall dimension w x l x h	355 x 380 x 415 mm		
Gewicht, netto	Net weight	16 kg		
Testglaseinsätze Typ 1 bis 4	tube racks typ 1 to 4	2 St.		
Bestell-Nr.	Order-No.	688.0032-S3	688.0033-S3	688.0034-S3

Technische Daten	Technical Data	B20-1	B20-2	B20-3
Arbeitstemperaturbereich max	Operating temperature range max		120°C	
Füllvolumen.	Bath capacity		20 lit.	
Badöffnung B x T / Badtiefe	Bath opening w x l / bath depth	300 x 370 / 150 mm		
Abmessungen B x T x H	Overall dimension w x l x h	355 x 558 x 375 mm		
Gewicht, netto	Net weight	17 kg		
Testglaseinsätze Typ 1 bis 4	tube racks typ 1 to 4	4 St.		
Bestell-Nr.	Order-No.	688.0035-S3	688.0036-S3	688.0037-S3

Technische Daten	Technical Data	B25-1	B25-2	B25-3
Arbeitstemperaturbereich max	Operating temperature range max		200°C	
Füllvolumen.	Bath capacity		25 lit.	
Badöffnung B x T / Badtiefe	Bath opening w x l / bath depth	300 x 370 / 200 mm		
Abmessungen B x T x H	Overall dimension w x l x h	355 x 558 x 415 mm		
Gewicht, netto	Net weight	19 kg		
Testglaseinsätze Typ 1 bis 4	tube racks typ 1 to 4	4 St.		
Bestell-Nr.	Order-No.	688.0038-S3	688.0039-S3	688.0040-S3

Zubehör	Accessoires:	B8	B12	B15	B20	B25
Entleerungshahn	drain cock	6026	6026	6026	6026	6026
Badbrücke	bath bridge	6303	6308	6308	6308	6308
variabl. Stellboden	var. Bottom		6297	6297	6298	6298
Bad-Deckel vorn	bath cover front	6214	6023	6023	6023	6023
Bad-Deckel hinten	bath cover back				6024	6024
Bad-Deckel einteilig	bath cover 1-piece				6025	6025

Testglaseinsätze / Test tube racks

aus Edelstahl für Bad B12-B25 / made of stainless steel for bath B12-B25

Typ / Type	Bohrungen / Bore holes	Eintauchtiefe / Depth [mm]	Best.Nr. / Order.No.
1	36 x ø17	100	6037
2	45 x ø13	70	6038
3	46 x ø17	100	6039
4	58 x ø13	70	6040

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Circulator bath with air-cooled refrigeration unit. Plug&Play Technology - microprocessor controlled, easy handling with encoder and digital display. Powerful force pump, moisted parts in stainless steel or high-resistant plastics. All models with 3 years warranty for electronic. 3 models available:

Polystat K6-1, with level protection (float switch) for continuous operation without personal assistance if using non-flammable liquids (FL). Maximum and minimum setpoint for additional safety. Max. temperature 200°C.

Polystat K6-2, similar to Polystat cc1, but with adjustable overtemperature protection (FL), programmer (5 steps), 3 fixtemperatures and easy ramping function. Max. temperature 200°C.

Polystat K6-3, with all functions of the Compatible Control Thermostats: Programmer (50 steps, divisible into 10 programs) temperature sequence controller for external thermoregulation and interfaces RS 232, RS 485 and analog (4...20mA) for bidirectional communication (FL). Serial printer output. Max. temperature 200°C.

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Accessoires and periphery: pump adapter for external thermoregulation*, adapter nom. dia 8 mm* / 12 mm*, dummy plugs* and sleeve nuts thread M16x1*, micro boltings, connection tubes, connection cable*. units for cells and glasses.

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Output data go for: room temperature 20°C

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The cooling baths are water baths with air-cooled refrigerating unit. Housing (K12 of plastic) and tank of stainless steel. The refrigerating unit works continuously.

*Together with the immersion thermostats can the cooling baths be used up to 200°C (K12 up to 120°C). The cooling unit can operate up to 50°C.

Technische Daten		Technical Data		Kältebad K12-S3
Arbeitstemperaturbereich		Operating temperature range		-20...(50) 120°C*
Kälteleistung	bei 0°C	Cooling capacity	at 0°C	0,18 kW
	bei -10°C		at -10°C	0,10 kW
Kältemittel		Refrigerent		R134a
Füllvolumen min.-max.		Bath capacity from - to		8 - 12 lit.
Badöffnung B x T / Badtiefe		Bath opening w x l / bath depth		295 x 320 / 150 mm
Abmessungen B x T x H		Overall dimension w x l x h		350 x 558 x 125 mm
Gewicht, netto		Net weight		20 kg
Netzanschluss		Power supply requirement		110-120 V ~ 50/60 Hz
Leistungsaufnahme/Absicherung		Power input / fuse		326 Watt 6 Amp.
Bestell-Nr.		Order-No.		653.0020-S3
ab Fert.Nr.		from Serial No.		52180 V1.3/02

Technische Daten		Technical Data		Kältebad K15-S3
Arbeitstemperaturbereich		Operating temperature range		-20...(50) 200°C*
Kälteleistung	bei 0°C	Cooling capacity	at 0°C	0,20 kW
	bei -10°C		at -10°C	0,12 kW
Kältemittel		Refrigerent		R134a
Füllvolumen min.-max.		Bath capacity from - to		11 - 15 l
Badöffnung B x T / Badtiefe		Bath opening w x l / bath depth		295 x 320 / 200 mm
Abmessungen B x T x H		Overall dimension w x l x h		355 x 558 x 240 mm
Gewicht, netto		Net weight		19 kg
Netzanschluss		Power supply requirement		110-120 V ~ 50/60 Hz
Leistungsaufnahme/Absicherung		Power input / fuse		326 Watt 6 Amp.
Bestell-Nr.		Order-No.		645.0021-S3
ab Fert.Nr.		from Serial.No.		46720 V1.2

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Bath- and circulation-thermostat. Housing, bath and all moisted parts in stainless steel. With cooling coil for water-cooling, pressure- and suction pump, level protection (float switch) and adjustable overtemperature protection.

PLUG & PLAY Technology: Microprocessor controller COMPATIBLE CONTROL with 3-year warranty. The controllers are interchangeable crosswise and can be calibrated. 3 models are available.

CC302-1, with level protection (float switch) and adjustable overtemperature protection for continuous operation without personal assistance for using flammable and non-flammable liquids (FL). Maximum and minimum setpoint for additional safety.

CC302-2, similar to CC302-1, but with programmer (5 steps), temperature sequence controller for external thermoregulation, 3 fixtemperatures and easy ramping function.

CC302-3, similar to CC302-2, but with all functions of the Compatible Control Thermostats: Programmer (50 steps, divisible into 10 programmes) and interfaces RS232, RS485 and analog (4...20mA) for bidirectional communication. Serial printer output.

Technische Daten		Technical Data	CC302-1	CC302-2	CC302-3
Arbeitstemperaturbereich		Operating temperature range		60...300 °C	
mit Wasserkühlung		with water cooling		20...300 °C	
mit Kühlgerät		with refrigeration chiller		-20...300 °C	
Temperaturkonstanz bei 70 °C		Temperature stability at 70 °C		0,02 K (DIN 58966)	
Temperatureinstellung		Temperature adjustment		digital	
Temperaturanzeige		Temperature indication		digital	
Absolutgenauigkeit		absolute accuracy	kalibrierbar / setup for calibration		
Temperaturfühler		Temperature sensor		Pt 100	
externer Programmeingang		external program input	---	---	4-20mA
Schreiber Ausgang		Recorder output	---	---	4-20mA
Sicherheitsklasse		Safety classification		FL	
Heizleistung		Heating capacity		2,0 kW	
Druckpumpe	Anschluß 12mm	Force pump (adapter nom 12 mm)		18 l/min max.	
	Anschluß 8 mm	(adapter nom 8 mm)			
Förderhöhe (Druck)		Pressure		0,5 bar max.	
Saugpumpe		Suction pump		15 l/min max.	
Saughöhe (Sog)		Suction		0,4 bar max.	
Pumpenanschluß		Pump connection		M16x1	
Füllvolumen max.		Bath capacity max.		8,5 lit.	
mit Verdrängereinsatz		with displacement rack		5,2 lit.	
Badöffnung B x T / Tiefe		Bath opening w x l / depth		180 x 120 / 155 mm	
Abmessungen B x T x H		Overall dimension w x l x h		230 x 365 x 405 mm	
Arbeitshöhe Bad		Heith of bath opening		190 mm	
Gewicht, netto		Net weight		18,2 kg	
Netzanschluß		Power supply requirement		230 V ~ 50 Hz	
Absicherung		fuse		10 A	
Bestell-Nr.		Order-No.	683.0001	683.0002	683.0003
ab Fert.Nr.		from serial no.		36180	V1.0

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Accessoires and periphery: Adapter nom. dia 12 mm*, dummy plugs* and sleeve nuts thread M16x1*, adapter nom. dia 8 mm, micro boltings, connection tubes, displacement rack for reducement of the bath volume and temperature sequence control for external thermoregulation, digital programmers, serial interface and current voltage converter. * standard equipment

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


(from Software V3.04)

Features of the set-key

The set key is for choosing the set point function or the menu function and, in general, for confirming an input. Changes and new settings have to be confirmed resp. memorized.

In case of a flashing display, the Set key serves as Enter key.

Enter =  when display is blinking

Call up the Menu functions

Press the set key and switch on the thermostat at the mains switch.
Wait is displayed. Release the set key. After approx. 5 Sec. the first function of the function menu Limit is displayed. You can leaf through the different functions of the function menu using the encoder. By pressing the set key, the chosen function is called up.

Setting Parameters

With the help of the  (encoder) the chosen function can be activated.