


Progress on RICH detector

Performance Estimation and Technical Issues

Rome/Sanità Group

11 April 2006 / 2nd Transversity Meeting

RICH schematic view




- Radiator: Liquid Freon (15 mm) at Room temperature
- Radiator Vessel: $1.9 \times 0.31 \times$ Neoceram (4 mm) / Quartz (5 mm)
- Proximity Gap: 10 cm
- Photon Detector: MWPC + 3 PAD PLANES (~ 11000 pads)
- Dry Methane at STP


+ the evaporation facility
(now at Stony Brook University)

Current RICH performance at $2.0 \pm 5\%$ GeV/c $\Rightarrow \pi/K$ rejection 1:1000

- Mean Number of Firing Pads (for π) : $N_{PAD,\pi} = 13$
- Cherenkov Angle resolution (sigma) : $\sigma_{\theta_{\pi}} = 5$ mrad
 $\rightarrow \Delta\theta_{\pi,K} = 6\sigma_{\theta_{\pi}}$



Pattern Recognition



Pretty good agreement between experimental data and MC:

$$\sigma_{\theta_{\pi}}^{EX} / \sigma_{\theta_{\pi}}^{MC} \sim 0.9$$

$$N_{PAD,\pi}^{EX} / N_{PAD,\pi}^{MC} \sim 0.9$$

$$\Delta\theta_{\pi,K}^{EX} / \Delta\theta_{\pi,K}^{MC} \sim 1.07$$

2.4 GeV/c Upgrade

Requirement (?)

Rule of thumb: π/K rejection factor at the level of 1:1000

$$\Rightarrow \Delta\theta_{\pi,K} \geq 5.5 \cdot \sigma_{\theta}$$

Criteria

- Re-use existing (spare) parts as much as possible
- Minimize costs and resources
- Guarantee the existing performance (switching back to the pre-upgrade conditions in few hours)

Constraints

- Maximum available electronics channels: 2×11520
- Maximum pad plane size: $65 \times 45 \text{ cm}^2$ (evaporation chamber)
- Maximum size along particle trajectory $\sim 1 \text{ m}$
- ???

... HRS Electron Detectors Layout


RICH installed at 2.6 m (from Bottom VDC) in HRS-Hadron

... HRS Electron Detectors Layout

Electron Arm

Detector

Package




... HRS Electron Detectors Layout


Electron Arm

Detector


Package



... HRS Electron Detectors Layout



... HRS Electron Detectors Layout



can be moved to ~ 1.7 m in HRS-Left

2.4 GeV/c Upgrade Options

Optimize the **gap length** and ...

- ① **NO CHANGES**: leave as it is
- ② **NEW FREON**: use radiator with lower refractive index (C_5F_{12} instead of C_6F_{14})
- ③ **2×PLANE**: double the photon detection surface (double the pad planes)
- ④ **PAD/2**: halve the pad size (double the pad numbers) on the same pad plane
- ⑤ **MIXED**: quite a few combinations of the above; e.g.:
 - ▶ **NEW RADIATOR** and **2×PLANE**
 - ▶ Reduce the pad size and increase the detection surface appropriately (up to the maximum available channels)

Options 3 and 4 are exclusive (max. electronics channel)

2.4 GeV/c Upgrade Options

Optimize the **gap length** and ...

- ① NO CHANGES: leave as it is
- ② **NEW FREON**: use radiator with lower refractive index (C_5F_{12} instead of C_6F_{14})
- ③ $2\times$ PLANE: double the photon detection surface (double the pad planes)
- ④ PAD/2: halve the pad size (double the pad numbers) on the same pad plane
- ⑤ MIXED: quite a few combinations of the above; e.g.:
 - ▶ NEW RADIATOR and $2\times$ PLANE
 - ▶ Reduce the pad size and increase the detection surface appropriately (up to the maximum available channels)

Options 3 and 4 are exclusive (max. electronics channel)

2.4 GeV/c Upgrade Options

Optimize the **gap length** and ...

- ① NO CHANGES: leave as it is
- ② NEW FREON: use radiator with lower refractive index (C_5F_{12} instead of C_6F_{14})
- ③ **2×PLANE**: double the photon detection surface (double the pad planes)
- ④ PAD/2: halve the pad size (double the pad numbers) on the same pad plane
- ⑤ MIXED: quite a few combinations of the above; e.g.:
 - ▶ NEW RADIATOR and 2×PLANE
 - ▶ Reduce the pad size and increase the detection surface appropriately (up to the maximum available channels)

Options 3 and 4 are exclusive (max. electronics channel)

2.4 GeV/c Upgrade Options

Optimize the **gap length** and ...

- ① NO CHANGES: leave as it is
- ② NEW FREON: use radiator with lower refractive index (C_5F_{12} instead of C_6F_{14})
- ③ 2×PLANE: double the photon detection surface (double the pad planes)
- ④ PAD/2: halve the pad size (double the pad numbers) on the same pad plane
- ⑤ MIXED: quite a few combinations of the above; e.g.:
 - ▶ NEW RADIATOR and 2×PLANE
 - ▶ Reduce the pad size and increase the detection surface appropriately (up to the maximum available channels)

Options 3 and 4 are exclusive (max. electronics channel)

2.4 GeV/c Upgrade Options



Optimize the **gap length** and ...

- ① **NO CHANGES**: leave as it is
- ② **NEW FREON**: use radiator with lower refractive index (C_5F_{12} instead of C_6F_{14})
- ③ **2×PLANE**: double the photon detection surface (double the pad planes)
- ④ **PAD/2**: halve the pad size (double the pad numbers) on the same pad plane
- ⑤ **MIXED**: quite a few combinations of the above; e.g.:
 - ▶ **NEW RADIATOR** and **2×PLANE**
 - ▶ Reduce the pad size and increase the detection surface appropriately (up to the maximum available channels)

Options 3 and 4 are exclusive (max. electronics channel)

Option 1 : NO CHANGES (Detector at 2.6 m from VDC)



Pattern Recognition



$$\Delta\theta_{\pi, K}^{MC} \sim 4.1 \cdot \sigma_\theta \Rightarrow \pi : K \sim 1 : 140 \text{ (95\% efficiency)}$$

Option 1 : NO CHANGES (Detector at 1.8 m from VDC)

Pattern Recognition





$$\Delta\theta_{\pi, K}^{MC} \sim 4.4 \cdot \sigma_\theta \Rightarrow \pi : K \sim 1 : 500 \text{ (95\% efficiency)}$$

Moving at 1.8 improve performance (better light collection)

Option 2 : NEW FREON (Detector at 1.8 m from VDC)

Pattern Recognition



$$\Delta\theta_{\pi, K}^{\text{MC}} \sim 5.5 \cdot \sigma_\theta \Rightarrow \pi : K < 1 : 1000 \text{ (95\% efficiency)}$$

Option 2 : NEW FREON: Technical/Thermal Issues

!! C5F12 boiling point @ 1atm [C] = 30 °C !!

Thermal Effects / Good News

- Very low neoceram (and quartz) expansion coefficient
- Epoxy elasticity will absorb dilatation effects (expected small)
- Mechanical stresses under control at $T = 21^{\circ}\text{C} \pm 20^{\circ}\text{C}$

Option 2 : NEW FREON: Technical/Thermal Issues

!! C5F12 boiling point @ 1atm [C] = 30 °C !!

Thermal Effects / Good News

- Very low neoceram (and quartz) expansion coefficient
- Epoxy elasticity will absorb dilatation effects (expected small)
- Mechanical stresses under control at $T = 21^{\circ}\text{C} \pm 20^{\circ}\text{C}$


Thermal Effects / Bad News

- Cavitation (?): hot spots in the liquid can cause localized high pressure ... but inclination (that is pressure variation) may compensate the decrease of vapor pressure with temperature
- Disadvantageous Thermal Conductivity: any cooling system **MUST NOT** touch the radiator vessel in order to avoid propagation of possible vibrations to the fragile vessel
- Turbulence and Instability: higher freon flux may generate turbulence and mechanical instability of the radiator vessel


NEW FREON: Cooling System

Lateral cooling rods along the vessel long-side is not effective.

→ Large cooling plate between the RICH entrance window and the neoceram vessel



NEW FREON: Cooling Stainless Steel Plate System




- Two foils of stainless steel 0.5 ÷ 1.0 mm thick
- Liquid cooler: 5 ÷ 10 mm thick
- One plane, size: $\sim 300 \times 1500 \text{ mm}^2$
- Weight: $< 10 \text{ Kg}$

Additional material on particle path: $\sim 2 \text{ g/cm}^2$ (now 1.3 g/cm^2)

Mechanically compatible with the existing structure (planarity to be verified)

French Company DATE (DEVELOPPMENT ET APPLICATION DES TECHNIQUES DE L'ENERGIE)

NEW FREON: Is a cooling system really needed ?




- Max freon speed w/o cooling (~ 30 min radiator vessel refill time / current system works at ~ 6 h) seems to be compatible with laminar flow in pipe ($R_c \sim 2200$) and vessel ($R_c > 10000$)
- Inlet freon temperature can be reduced to $10^\circ C$
- Coarse cooling (down to $30^\circ C$) of the external radiator will help)
- Finite element analysis for mechanical stability

Open Problems:


- Hot Spot
- Max Freon Speed

Option 3 : 2×PLANE (Detector at 1.8 m from VDC)

Pattern Recognition



rz/c6f14_2.28p2.52_3.5_0.80x0.84_1.80/ss_out.txt





$$\Delta\theta_{\pi,K}^{MC} \sim 6.0 \cdot \sigma_\theta \Rightarrow \pi : K \ll 1 : 1000 \text{ (95\% efficiency)}$$

2×PLANE

- NO relevant technical issues
- Design is straightforward
- Re-use existing pad planes
- Use two acquisition systems (old and new one, both available)
- Some more contamination from delta electrons due to the middle break (~ 3.5 cm). (however electron MIP - kaon trajectory greater than ~ 6 cm!!)
- ~ 15 cm longer detector
- Major Impacts:
 - ▶ New MWPC (largest cost)
 - ▶ Partially new mechanical supports for the front-end cards
 - ▶ Modify (increase) the glove box used in the evaporation of the CsI on the pad planes
 - ▶ Longer evaporation procedure
- Costs: quotation requested



Option 4 : PAD/2 (Detector at 1.8 m from VDC)



$$\Delta\theta_{\pi, K}^{MC} \sim 5.0 \cdot \sigma_\theta \Rightarrow \pi : K \leq 1 : 1000 \text{ (95\% efficiency)}$$

High front-end electronics compactness is demanding

Option 5 : NEW FREON & 2×PLANE (Detector at 1.8 m from VDC)




$$\Delta\theta_{\pi, K}^{\text{MC}} \sim 8.0 \cdot \sigma_\theta \Rightarrow \pi : K \lll 1 : 1000 \text{ (99\% efficiency)}$$

Phase Space


HRS Phase Space

Trajectory Phase Space at $z=1.8$ m



RICH PE/MIP Distribution

rz/c5f12_2.28p2.52_3.5_0.80x0.84_1.80/proxie_13.00



Pro's and Con's

Option	$\Delta\theta$	Pro	Con
1	4.4	Trivial Implementation Negligible cost and effort	Poor performance
2	5.5	Good performance	Open issues Involve the fragile vessel Slightly complex design Add material in path (?)
3	6.0	Very Good performance Straight Design	More space required Two readout systems Expensive (?)
4	5.0	Moderate performance	Electronics compactness Two readout systems
5	8.0	Superior performance	See 2 and 3

... than ?

- NO CHANGES and detector at 1.8 m: π/K rejection $\sim 1 : 500$ (not so bad)
- NEW FREON requires more analysis for cooling
- 2×PLANE better solution, but cost may be high (expecting quotation today!)
- PAD/2, complexity/benefit not convenient
- NEW FREON & 2×PLANE: superior performance, but ...

coming ...

- In the next couple of weeks expected results from the finite element analysis by the Engineer that designed the radiator vessel
- Careful evaluation of the 2×PLANE to reduce costs