



PF-48CD-48x10


Connection scheme

FTSH con- nec- tor	Detector	FTSH		FTSH		FTSH		FTSH	
		FTSH	FTSH	FTSH	FTSH	FTSH	FTSH	FTSH	FTSH
1	mass	39	73	81	153	123	42	165	122
2	Floating	40	75	82	155	124	44	166	124
3	1	41	77	83	157	125	46	167	126
4	3	42	79	84	159	126	48	168	128
5	5	43	81	85	161	127	50	169	130
6	7	44	83	86	163	128	52	170	132
7	9	45	85	87	165	129	54	171	134
8	11	46	87	88	167	130	56	172	136
9	13	47	89	89	169	131	58	173	138
10	15	48	91	90	171	132	60	174	140
11	17	49	93	91	173	133	62	175	142
12	19	50	95	92	175	134	64	176	144
13	21	51	mass	93	177	135	66	177	146
14	23	52	mass	94	179	136	68	178	148
15	25	53	97	95	181	137	70	179	150
16	27	54	99	96	183	138	72	180	152
17	29	55	101	97	185	139	74	181	154
18	31	56	103	98	187	140	76	182	156
19	33	57	105	99	189	141	78	183	158
20	35	58	107	100	191	142	80	184	160
21	37	59	109	101	mass	143	82	185	162
22	39	60	111	102	GR	144	84	186	164
23	41	61	113	103	2	145	86	187	166
24	43	62	115	104	4	146	88	188	168
25	45	63	117	105	6	147	90	189	170
26	47	64	119	106	8	148	92	190	172
27	49	65	121	107	10	149	94	191	174
28	51	66	123	108	12	150	96	192	176
29	53	67	125	109	14	151	mass	193	178
30	55	68	127	110	16	152	mass	194	180
31	57	69	129	111	18	153	98	195	182
32	59	70	131	112	20	154	100	196	184
33	61	71	133	113	22	155	102	197	186
34	63	72	135	114	24	156	104	198	188
35	65	73	137	115	26	157	106	199	190
36	67	74	139	116	28	158	108	200	192
37	69	75	141	117	30	159	110		
38	71	76	143	118	32	160	112		
		77	145	119	34	161	114		
		78	147	120	36	162	116		
		79	149	121	38	163	118		
		80	151	122	40	164	120		



CANBERRA
PF-192CD-48*10





PF-48CD-48*10		Name: CERAMIC BOARD	
Date: 22/12/06		Doc: STTPR497B	
MAT: Al ₂ O ₃		Auth.: MAT	
SCALE: 2/1		Order code: SOR6339	
TOL: +/-0.1		Internal code: HW	